


NANOVA
AS-Q-CLIP[®]

***Non-absorbable Polymer
Ligation Clip System***



Polymer Clip with Patented Tissue-Locking Mechanism

NANOVA[™]

Safe and Secure Closure

BETTER BY DESIGN

Proprietary anti-migration features, or teeth, provide strong migration resistance reducing the likelihood of post-op clip migration, thus increasing surgeon's peace of mind.



ZERO LEARNING CURVE

Surgeons can enjoy the peace of mind that comes from the superior performance of Vas-Q-Clips® right away. There are no procedural changes when switching from other non-absorbable polymer ligation clips.



ML, L, AND XL

Enjoy the quality and performance of Vas-Q-Clip® clips in a variety of procedures with our 3 different sizes.





CLIP SECURITY


Safe and secure clip closure during application provides tactile feedback, assuring the surgeon that the Vas-Q-Clip® has encompassed the ligated vessel.



VESSEL SIZE RANGE



THE VAS-Q-CLIP® DIFFERENCE



DESIGN

The Vas-Q-Clip® is a direct result of Nanova Biomaterials Inc.'s mission. In the design phase, Nanova Biomaterials Inc. identified migration resistance and occlusion as key safety factors. When considering clip migration and occlusion, the Vas-Q-Clip® withstands a high force and pressure respectively. Patient Safety is our top priority.

The Vas-Q-Clip® is designed and manufactured in the USA

VAS-Q-CLIP® LIGATION CLIPS



Catalog Reference	Clip Size	Vessel Range	Clips/Cartridge	Cartridge/Box
2231	Medium-Large	3-10mm	6	14
2232	Large	5-13mm	6	14
2233	Extra Large	7-16mm	6	14

ENDOSCOPIC CLIP APPLIERS



Catalog Reference	Size	Diameter	Length	Jaw Angle
1161-ML	Medium-Large	5mm	33cm	Straight
1162-L	Large	10mm	33cm	Curved
1163-XL	Extra Large	10mm	33cm	Curved
1461-ML	Medium-Large	5mm	45cm	Straight
1462-L	Large	10mm	45cm	Curved
1463-XL	Extra Large	10mm	45cm	Curved

OPEN CLIP APPLIERS



Catalog Reference	Size	Length	Jaw Angle
1261-ML	Medium-Large	27cm	Curved
1262-L	Large	27cm	Curved
1263-XL	Extra Large	27cm	Curved

STERILIZATION/STORAGE TRAYS



ID	Tray Capacity	Appliers Included	Appliers Working Length	Quantity
1879-Tray	3 Appliers	None	N/A	1
1879-116-Tray-Set	3 Appliers	ML, L, XL	32.5cm (12 3/4")	1

TO LEARN MORE CONTACT YOUR NANOVA REPRESENTATIVE